

Product Datasheet

Goat FSHB ELISA Kit orb403697

Description	Goat follicle-stimulating hormone,FSH ELISA Kit
Reactivity	Goat
Range	20 mIU/ml-800 mIU/ml
Sensitivity	20mIU/ml
Tested Applications	ELISA
Target	FSHB
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm