
Product Datasheet

Goat GnRH ELISA Kit (orb403700)

Description	Goat gonadotropin-releasing hormone, GnRH ELISA Kit
Reactivity	Goat
Range	2.5 pg/mL-500 pg/mL
Sensitivity	1.25 pg/mL
Tested Applications	ELISA
Target	GnRH
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	D3YJ02
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm