

## Product Datasheet

### Goat GnRH ELISA Kit orb403700

<b>Description</b>	Goat gonadotropin-releasing hormone,GnRH ELISA Kit
<b>Reactivity</b>	Goat
<b>Range</b>	2.5 pg/ml-500 pg/ml
<b>Sensitivity</b>	1.25pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	GnRH
<b>Note</b>	For research use only.
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	<ul style="list-style-type: none"><li>• <b>Sample Type:</b> serum, plasma, tissue homogenates</li><li>• <b>Assay Time:</b> 1-5 hours</li><li>• <b>Sample Volume:</b> 50-100 <math>\mu</math>l</li><li>• <b>Detection Wavelength:</b> 450 nm</li></ul>