

Product Datasheet

Goat LEP ELISA Kit orb403721

Description	Goat Leptin,LEP ELISA kit
Reactivity	Goat
Range	0.625 ng/ml-40 ng/ml
Sensitivity	0.156 ng/ml
Tested Applications	ELISA
Target	LEP
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm