

Product Datasheet

Goat SAA ELISA Kit orb403726

Description	Goat Serum amyloid A(SAA) ELISA Kit
Reactivity	Goat
Range	100 ng/ml-8000 ng/ml
Sensitivity	125 ng/ml
Tested Applications	ELISA
Target	SAA
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm