

Product Datasheet

Goat IGF1 ELISA Kit orb403728

Description	Goat insulin-like growth factor 1(IGF-1)ELISA Kit
Reactivity	Goat
Range	80 ng/ml -1000 ng/ml
Sensitivity	20 ng/ml
Tested Applications	ELISA
Target	IGF1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm