

## Product Datasheet

### Goat ADP ELISA Kit (orb403731)

<b>Description</b>	Goat adiponectin(ADP)ELISA Kit...
<b>Reactivity</b>	Goat
<b>Range</b>	0.312 ng/mL-20 ng/mL
<b>Sensitivity</b>	0.078 ng/mL
<b>Target</b>	ADP
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>U6A087</b>
<b>Expiration Date</b>	Please enquire.