

Product Datasheet

Guinea Pig INS ELISA Kit orb403742

Description	Guinea pig insulin (INS) ELISA kit
Reactivity	Guinea pig
Range	12 µIU/ml-210 µIU/ml
Sensitivity	7.5 µIU/ml
Tested Applications	ELISA
Target	INS
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm