

Product Datasheet

Guinea Pig IL10 ELISA Kit orb403750

Description	Guinea pig interleukin 10,IL-10 ELISA Kit
Reactivity	Guinea pig
Range	31.25 pg/ml-2000 pg/ml
Sensitivity	7.81pg/ml
Tested Applications	ELISA
Target	IL10
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm