

Product Datasheet

Guinea Pig IL1B ELISA Kit orb403755

Description	Guinea pig Interleukin-1 beta(IL1B) ELISA kit
Reactivity	Guinea pig
Range	15.6 pg/mL-1000 pg/mL
Sensitivity	3.9 pg/mL
Tested Applications	ELISA
Target	IL1B
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9WVG1
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm