

Product Datasheet

Guinea Pig IL5 ELISA Kit orb403760

Description	Guinea pig interleukin 5,IL-5 ELISA Kit
Reactivity	Guinea pig
Range	31.25 pg/mL-2000 pg/mL
Sensitivity	7.81 pg/mL
Tested Applications	ELISA
Target	IL5
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O08987
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma and tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm