

## Product Datasheet

### Guinea Pig MMP-9 ELISA Kit (orb403769)

<b>Description</b>	guinea pig matrix metalloproteinase 9,MMP-9 ELISA Kit...
<b>Reactivity</b>	Guinea pig
<b>Range</b>	3.12 ng/mL-200 ng/ mL
<b>Sensitivity</b>	0.78 ng/mL
<b>Target</b>	MMP-9
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>H0VAD9</b>
<b>Expiration Date</b>	Please enquire.