

Product Datasheet

Guinea Pig tri-iodothyronine ELISA Kit orb403779

Description	Guinea pig tri-iodothyronine (T3) ELISA kit
Reactivity	Guinea pig
Range	0.5ng/ml-8ng/ml
Sensitivity	0.5 ng/ml
Tested Applications	ELISA
Target	tri-iodothyronine
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm