

Product Datasheet

Guinea Pig TNF ELISA Kit orb403782

Description	Guinea pig Tumor necrosis factor α , TNF- α ELISA Kit
Reactivity	Guinea pig
Range	3.12 pg/mL-200 pg/mL
Sensitivity	0.78 pg/mL
Tested Applications	ELISA
Target	TNF
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P51435
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 μ l • Detection Wavelength: 450 nm