

## Product Datasheet

### Hamster CRP ELISA Kit orb403786

<b>Description</b>	Hamster C-reactive protein(CRP) ELISA Kit
<b>Reactivity</b>	Hamster
<b>Range</b>	12.5 ng/ml-800 ng/ml
<b>Sensitivity</b>	3.12ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	CRP
<b>Note</b>	For research use only.
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	<ul style="list-style-type: none"><li>• <b>Sample Type:</b> serum, plasma, tissue homogenates</li><li>• <b>Assay Time:</b> 1-5 hours</li><li>• <b>Sample Volume:</b> 50-100 <math>\mu</math>l</li><li>• <b>Detection Wavelength:</b> 450 nm</li></ul>