

## Product Datasheet

### Hamster IFN gamma ELISA Kit orb403788

<b>Description</b>	Hamster interferon, gamma (IFNG) ELISA kit
<b>Reactivity</b>	Hamster
<b>Range</b>	1.87 pg/mL-120 pg/mL
<b>Sensitivity</b>	0.46 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	IFNG
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenate
<b>Uniprot ID</b>	<b>O35497</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, tissue homogenate • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm