

## Product Datasheet

### GNPTAB antibody {HRP} orb157234-HRP

<b>Description</b>	Rabbit polyclonal antibody against GNPTAB conjugated to HRP.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Predicted Reactivity</b>	Porcine
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ICC, IF
<b>Immunogen</b>	KLH conjugated synthetic peptide derived from the region aa11-90 of human GNPTAB
<b>Target</b>	GNPTAB {HRP}
<b>Form/Appearance</b>	Liquid: PBS and 2mM Sodium Azide
<b>Concentration</b>	100 µg in 500 µl (0.20 mg/ml)
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purified by Protein A
<b>MW</b>	105 kDa
<b>Uniprot ID</b>	<b>Q3T906</b> (Human)
<b>Entrez</b>	<b>79158</b>
<b>Dilution Range</b>	To be determined by the end user
<b>Expiration Date</b>	12 months from date of receipt.