

## Product Datasheet

### Rapid APC Antibody Labeling Kit (orb867211)

<b>Description</b>	Rapid APC Antibody Labeling Kit...
<b>Conjugation</b>	APC
<b>Tested Applications</b>	ELISA, FC, IHC-P, WB
<b>Hazard Information</b>	Non-Toxic
<b>Storage</b>	Store at (2-8°C), Minimize light exposure
<b>Note</b>	For research use only.
<b>Dilution Range</b>	To Be determined by end-user
<b>Expiration Date</b>	12 months from date of receipt.