

Product Datasheet

Horse INS ELISA Kit orb403816

Description	Horse insulin (INS) ELISA kit
Reactivity	Equine
Range	1.56 nIU/ml-100 nIU/ml
Sensitivity	0.39 nIU/ml.
Tested Applications	ELISA
Target	INS
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm