

Product Datasheet

Horse IFN gamma ELISA Kit [Out of Stock] orb403819

Description	Horse Interferon gamma(IFNG) ELISA kit
Reactivity	Equine
Range	6.25 pg/ml-400 pg/ml
Sensitivity	1.56pg/ml
Tested Applications	ELISA
Target	IFN gamma
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm