

Product Datasheet

Horse LEP ELISA Kit orb403834

Description	Horse Leptin(LEP) ELISA kit
Reactivity	Equine
Range	3.12 pg/ml-200 pg/ml
Sensitivity	0.78 pg/ml
Tested Applications	ELISA
Target	LEP
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm