

Product Datasheet

Horse Prolactin ELISA Kit orb403854

Description	Horse Prolactin(PRL) ELISA kit
Reactivity	Equine
Range	2.5 ng/ml-1000 ng/ml
Sensitivity	2.5 ng/ml
Tested Applications	ELISA
Target	Prolactin
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm