

Product Datasheet

Horse HIF1A ELISA Kit orb403859

Description	Horse hypoxia inducible factor 1, alpha subunit (basic helix-loop-helix transcription factor) (HIF1A) ELISA kit
Reactivity	Equine
Range	31.25 pg/ml-2000 pg/ml
Sensitivity	7.81 pg/ml
Tested Applications	ELISA
Target	HIF1A
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm