

Product Datasheet

Horse TGF beta 1 ELISA Kit orb403868

Description	Horse Transforming growth factor beta-1(TGFB1) ELISA kit
Reactivity	Equine
Range	3.12 pg/mL-200 pg/mL
Sensitivity	1.56 pg/mL
Tested Applications	ELISA
Target	TGFB1
Note	For research use only.
Sample Types	serum, plasma
Uniprot ID	O19011
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm