

Product Datasheet

Horse PMSG ELISA Kit orb403869

Description	Pregnant mare serum gonadotrophin (PMSG)ELISA Kit
Reactivity	Equine
Range	5 mIU/mL-2000 mIU/mL
Sensitivity	5 mIU/mL
Tested Applications	ELISA
Target	PMSG
Note	For research use only.
Sample Types	serum
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm