

Product Datasheet

Human PSMD4 ELISA Kit orb403879

Description	Human 26S proteasome non-ATPase regulatory subunit 4(PSMD4) ELISA kit
Reactivity	Human
Range	34.4 pg/ml-2200 pg/ml
Sensitivity	8.6 pg/ml
Tested Applications	ELISA
Target	PSMD4
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm