

## Product Datasheet

### Human ADAMTS14 ELISA Kit orb403894

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human A disintegrin and metalloproteinase with thrombospondin motifs 14(ADAMTS14) ELISA kit   |
| <b>Reactivity</b>          | Human   |
| <b>Tested Applications</b> | ELISA   |
| <b>Target</b>              | ADAMTS14  |
| <b>Note</b>                | For research use only.  |
| <b>Expiration Date</b>     | Please enquire.   |
| <b>Application Notes</b>   | <ul style="list-style-type: none"><li>• <b>Assay Time:</b> 1-5 hours</li><li>• <b>Sample Volume:</b> 50-100 <math>\mu</math>l</li><li>• <b>Detection Wavelength:</b> 450 nm</li></ul> |