

Product Datasheet

Human AOA_H ELISA Kit orb403924

Description	Human Acyloxyacyl hydrolase(AOA _H) ELISA kit
Reactivity	Human
Range	0.156 ng/ml-10 ng/ml
Sensitivity	0.039 ng/ml
Tested Applications	ELISA
Target	AOA _H
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm