
Product Datasheet

Human AFM ELISA Kit (orb403939)

Description	Human Afamin(AFM) ELISA kit
Reactivity	Human
Range	3.9 ng/mL-250 ng/mL
Sensitivity	0.97 ng/mL
Tested Applications	ELISA
Target	AFM
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P43652
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm