

Product Datasheet

Human ORM1 ELISA Kit orb403951

Description	Human Alpha-1-acid glycoprotein 1(ORM1) ELISA kit
Reactivity	Human
Range	31.25 ng/mL-2000 ng/mL
Sensitivity	7.8 ng/mL
Tested Applications	ELISA
Target	ORM1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P02763
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm