

Product Datasheet

Human ORM1 ELISA Kit orb403951

Description	Human Alpha-1-acid glycoprotein 1(ORM1) ELISA kit
Reactivity	Human
Range	31.25 ng/ml-2000 ng/ml
Sensitivity	7.8 ng/ml
Tested Applications	ELISA
Target	ORM1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm