

Product Datasheet

Human DEFA7 ELISA Kit orb403956

Description	Human Alpha-defensin 7(DEFA7) ELISA kit
Reactivity	Human
Range	0.625 ng/ml-40 ng/ml
Sensitivity	0.156 ng/ml.
Tested Applications	ELISA
Target	DEFA7
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm