

Product Datasheet

Human SCYE1 ELISA Kit (orb403963)

Description	Human Aminoacyl tRNA synthase complex-interacting multifunctional protein 1(SCYE1) ELISA kit
Reactivity	Human
Range	23.5 pg/mL-1500 pg/mL
Sensitivity	5.8 pg/mL
Tested Applications	ELISA
Target	SCYE1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q12904
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm