

## Product Datasheet

### Human KAL1 ELISA Kit (orb403978)

<b>Description</b>	Human Anosmin-1(KAL1) ELISA kit...
<b>Reactivity</b>	Human
<b>Range</b>	125 pg/mL-8000 pg/mL
<b>Sensitivity</b>	31.25 pg/mL
<b>Target</b>	KAL1
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">P23352</a>
<b>Expiration Date</b>	Please enquire.