

Product Datasheet

Human ASS1 ELISA Kit orb403994

Description	Human Argininosuccinate synthase(ASS1) ELISA kit
Reactivity	Human
Range	7.8 pg/ml-500 pg/ml
Sensitivity	1.95 pg/ml
Tested Applications	ELISA
Target	ASS1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm