

Product Datasheet

Human ASPN ELISA Kit orb404000

Description	Human Asporin(ASPN) ELISA kit
Reactivity	Human
Range	34.37 pg/ml-2200 pg/ml
Sensitivity	8.59 pg/ml
Tested Applications	ELISA
Target	ASPN
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm