

Product Datasheet

Human BTLA ELISA Kit orb404009

Description	Human B- and T-lymphocyte attenuator(BTLA) ELISA kit
Reactivity	Human
Range	0.45 ng/ml-30 ng/ml
Sensitivity	0.11 ng/ml
Tested Applications	ELISA
Target	BTLA
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm