

## Product Datasheet

### Human BCL2L11 ELISA Kit orb404018

<b>Description</b>	Human Bcl-2-like protein 11(BCL2L11) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	23.44 pg/mL-1500 pg/mL
<b>Sensitivity</b>	5.86 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	BCL2L11
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<b>O43521</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm