

Product Datasheet

Human ADRBK1 ELISA Kit orb404022

Description	Human Beta-adrenergic receptor kinase 1(ADRBK1) ELISA kit
Reactivity	Human
Range	18.75 pg/ml-1200 pg/ml
Sensitivity	4.68 pg/ml
Tested Applications	ELISA
Target	ADRBK1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm