

Product Datasheet

Human BMP8A ELISA Kit orb404046

Description	Human Bone morphogenetic protein 8A(BMP8A) ELISA kit
Reactivity	Human
Range	37.5 pg/mL-2400 pg/mL
Sensitivity	9.38 pg/mL
Tested Applications	ELISA
Target	BMP8A
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q7Z5Y6
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, cell lysates. • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm