

Product Datasheet

Human BMP8B ELISA Kit orb404047

| | |
|----------------------------|---|
| Description | Human Bone morphogenetic protein 8B(BMP8B) ELISA kit |
| Reactivity | Human |
| Range | 9.37 pg/ml-600 pg/ml |
| Sensitivity | 2.34 pg/ml |
| Tested Applications | ELISA |
| Target | BMP8B |
| Note | For research use only. |
| Expiration Date | Please enquire. |
| Application Notes | <ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm |