

Product Datasheet

Human CPE ELISA Kit orb404082

Description	Human Carboxypeptidase E(CPE) ELISA kit
Reactivity	Human
Range	39.06 pg/ml-2500 pg/ml
Sensitivity	9.76 pg/ml.
Tested Applications	ELISA
Target	CPE
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm