
Product Datasheet

Human CPE ELISA Kit (orb404082)

Description	Human Carboxypeptidase E(CPE) ELISA kit
Reactivity	Human
Range	39.06 pg/mL-2500 pg/mL
Sensitivity	9.76 pg/mL
Tested Applications	ELISA
Target	CPE
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P16870
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm