
Product Datasheet

Human CPN2 ELISA Kit (orb404084)

Description	Human Carboxypeptidase N subunit 2(CPN2) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Sensitivity	0.039 ng/mL
Tested Applications	ELISA
Target	CPN2
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P22792
Expiration Date	Please enquire.
Application Notes	• Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm