

Product Datasheet

Human CRTAC1 ELISA Kit orb404095

Description	Human Cartilage acidic protein 1(CRTAC1) ELISA kit
Reactivity	Human
Range	4.68 ng/ml-300 ng/ml
Sensitivity	1.17 ng/ml
Tested Applications	ELISA
Target	CRTAC1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm