
Product Datasheet

Human CD160 ELISA Kit (orb404113)

Description	Human CD160 antigen(CD160) ELISA kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Sensitivity	0.156 ng/mL
Tested Applications	ELISA
Target	CD160
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O95971
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm