

Product Datasheet

Human CD2AP ELISA Kit orb404119

Description	Human CD2-associated protein(CD2AP) ELISA kit
Reactivity	Human
Range	18.75 pg/mL-1200 pg/mL
Sensitivity	4.7 pg/mL
Tested Applications	ELISA
Target	CD2AP
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9Y5K6
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm