

Product Datasheet

Human CENPK ELISA Kit orb404136

Description	Human Centromere protein K(CENPK) ELISA kit
Reactivity	Human
Range	23.5 pg/ml-1500 pg/ml
Sensitivity	5.8 pg/ml
Tested Applications	ELISA
Target	CENPK
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm