

Product Datasheet

Human CLDN3 ELISA Kit orb404155

Description	Human Claudin-3(CLDN3) ELISA kit
Reactivity	Human
Range	25 pg/ml-1600 pg/ml
Sensitivity	6.25 pg/ml
Tested Applications	ELISA
Target	CLDN3
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm