

Product Datasheet

Human C-MAPT/C-TAU ELISA Kit orb404157

Description	Human cleaved microtubule-associated protein tau (C-MAPT/C-TAU) ELISA kit
Reactivity	Human
Range	62.5 pg/ml-4000 pg/ml
Sensitivity	15.6 pg/ml.
Tested Applications	ELISA
Target	C-MAPT/C-TAU
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cerebrospinal fluid (CSF), tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450nm